Optimal V-Notched Shear Testing of Composite Laminae and Laminates

Dan Adams

Vice President
Wyoming Test Fixtures, Inc.
Salt Lake City, UT, USA
and
Emeritus Professor of Mechanical Engineering
University of Utah

Shear Strength Workshop March 7, 2022

Salt Lake City, UT, USA



Optimal Shear Test of Composites: What attributes come to mind?

Within the "test section":

- Uniform state of shear stress
- "Pure shear" no other stress components
- Higher magnitude of shear stress than other regions of the specimen
- Shear failure produced in test section



Comparison of V-Notched Shear Tests: Test Fixtures and Specimens



V-Notched Beam "Iosipescu" Shear ASTM D5379





V-Notched Rail Shear ASTM D7078







Combined Loading Shear



Iosipescu (V-Notched Beam) Shear Test: ASTM D5379

- 76 mm long x 19 mm wide specimen
- Opposing 90° machined V-notches
- Asymmetrical four-point flexure loading
- In-plane and interlaminar shear testing









V-Notched Rail Shear Test: ASTM D7078

- Standardized in 2005 by ASTM
- 76 mm x 56 mm notched specimen
- Same notch configuration as Iosipescu specimen
- Increased gage section (compared to Iosipescu)
- Face-loading allows testing of higher shear strength laminates









Combined Loading Shear Test (In Process of ASTM Standardization)


- Retains face loading of current D7078 V-notched rail shear test
- Introduce edge loading similar to D5379 losipescu shear test
- Specimen length increased to 127 mm
- Allow shear testing of thicker, higher shear strength composite laminates







Shear Stress Distribution in Test Section Area: V-Notched Shear Tests




Results from finite element simulation, AS4/3501-6 carbon/epoxy



Comparison of Loading Capabilities: V-Notched Rail Shear Tests

Shear load capability increased significantly using Combined Loading Shear (CLS) test method



Shear Failure of Quasi-Isotropic Combined Loading Shear (CLS) Test Specimen



IM7/8552 carbon/epoxy, ~10 mm specimen thickness

