



Join the Freezer Challenge!

Tips for efficient use of fridges and freezers:

Keep the door shut as much as possible.

Be clutter free; do not stack boxes against or on top of the units.

Report damaged seals immediately to the cold storage manager.

Use racking or equivalent storage solutions.

Keep current and visible inventory of items.

Unless specifically required, ULT freezers should be set at -70°C.

Cold Storage Maintenance

The Cold Storage Manager for this lab is:

.....

This unit will undergo maintenance every 6 months.

The preventative maintenance will include defrosting, inventory check, cleaning the filters, ensuring seals are intact, etc. If you have any questions or concerns, contact the Cold Storage Manager.



Cold Storage Maintenance Log

Date	Name

